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
User manual | Manuel d'utilisation | Istruzioni per l'uso | Manual de instrucciones

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SBS-LW-3000 | 2000 | 200

www.steinbergsystems.de

Inhalt

Contenu | Content | Contenuto | Contenido | Treść

Deutsch	4
English	8
Français	13
Italiano	18
Español	23
Polski	28

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
www.steinbergsystems.de

Summary

1. General Information
2. Warnings and notes
3. Package Contents
4. Technical specifications
5. Key Functions
6. Information service

1. General Information

The series of precision balances SBS-LW meets the latest technology of reliable processors with high-quality programs. By using the most current series of microprocessors AT and through integration of technologies for high speed and high precision, complemented by compensation circuits, these scales very quickly reach a stable weight mode. High anti-interference is also given. This series of precision balances are used very diversely and are characterized by their robust construction, particularly suited for use in industry. Moreover, it is very suitable for use in the commerce and in the laboratory. Our models are not calibrated, however. Therefore, they can only be used as instruments of control. The accuracy is still guaranteed by a calibration on the gravitational force average in Europe, allowing a low degree of tolerance.

2. Warnings and notes

The electronic weighing is a method of precise measurement. Security must remain at the forefront, please follow the instructions in this manual to increase the life of your device and ensure good long term results. Use the scale only under conditions for which it is conceived. Set it free from drafts, of varying temperatures and unstable environments, or electromagnetic.

3. Package Contents

Check by opening the package that all documents and accessories are included:


- Precision Balance
- Weighing Platform
- Power Cable
- Instructions

Please carefully read the instructions!

4. Technical specifications

Model	SBS-LW-3000 / 2000 / 200
Article	Precision Scale
Measuring range	3000/0,1g 2000/0,01g 200/0,001g
Dead weight	1,7kg
Functions	Recording of the automatic calibration Zeroing Tare Counting Error display
Environmental Conditions	5 ... +35 ° C / 10 ... 50 ° C (by use/ in stock) 50% / <- 85% humidity

5. Key Functions

WARNING SIGNALS AND STATUS

Overload (-----)

If overloaded, the device displays ----- (overload)

CAUTION

Discharge the balance immediately!

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By posting error I, the balance alerts you that the weight is mobile and unstable, and leads to measurement errors.

Instability (Error2)

By posting error2, the balance alerts you the weight is still unstable, which leads to errors of measurement.

On / Off (I / O)

With the main switch I / O on the back of the device, you can turn the scale on. A self-test begins. It begins with a screen test, which should display correctly 8.8.8.8.8.8.

then the maximum weight allowed by the balance (3000 g/2000g/200g).

After testing, the scale shows 0 and is ready for use.

To turn the scale off, press the button I / O again

Tare (TAR)

Place an empty container on the scale. Wait 3 sec. Until the scale displays the weight of the container. Press (TAR / CAL) to begin weighing ignoring the weight of the container.

CAUTION

The maximum load is reduced by the weight of the empty container!

For example:

Container b = 50g

Max. Normal l = 200g.

Max. remaining L = l - b = 150g.

Place the weight slowly on the scale and remove it immediately in case of overload (----- displayed)

Counting mode (COU)

In counting mode, you can determine the number of units put on the scale according to their weight. Place a sample on the scale. Switch to counting by pressing (COU / AFF).

With the (TAR / CAL) button, you can choose the number corresponding to the quantity of pieces placed (10 to 500 in increments of 10,20,50,100,250,500)

Confirm by pressing (COU / AFF). You can now weigh in counting mode.

calibration (CAL)

There are two ways to calibrate the scale.

The second way is more accurate, but requires the use of multiple standards.

1. Press the button (CAL) and put the weight that flashes on the scale. The display stabilizes, you can remove the standard. The scale is calibrated
2. Press and hold (COU / AFF) until the scale displays CAL. Place the weight corresponding to the flashing value on the scale.

Standards requirements for calibration

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
This Declaration of Conformity is suitable to the European Standard EN 45014 "General criteria for supplier's declaration of conformity". The basis for the criteria has been found in international documentation, particularly in: ISO/IEC Guide 22, 1982, "Informations on manufacturer's declaration of conformity with standards or other technical specifications".

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"Information on manufacturer's declaration of conformity with standards or other technical specifications".

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expondo GmbH